


Artificial Intelligence used to enhance decision making during colorectal cancer surgery for first time

17 December 2020



NIR intensities were extracted by tracking ROIs within each video, focusing on the initial ICG wash-in period of 100–300 s. Credit: UCD Research and Innovation

A clinical research study published in the *British Journal of Surgery* shows that fluorescence guidance (where a fluorescent chemical compound re-emits light with a long wavelength) can enable a colorectal surgeon to assess cancer tissues visually and with more specificity in real-time during surgery, by using near-infrared (NIR) light from an administered fluorophore in conjunction with artificial intelligence (AI) methods.

In this study, supported by the Disruptive Technologies and Innovation Fund 2018, videos from 24 patient (11 with cancer) surgeries were studied. Numerous ROIs (Regions of Interest) from each area of abnormality were selected for analysis from each video. NIR intensities were extracted by tracking the ROIs within each video, focusing on the initial wash-in period. The data set used for analysis included 435 ROI profiles each with 12 perfusion-characterizing features with balanced outcomes. At patient level, the system correctly diagnosed 19 of 20 cancers (95%).

Prof Ronan Cahill, professor of [surgery](#) at University College Dublin (UCD) and the Mater Misericordiae University Hospital (MMUH), said "Surgery has the substantial role to play in the therapy of over two-thirds of all cancers and key surgical decisions are traditionally made by human visual judgements, which assume a static biological FOV (Field of View) during the time frame of the observation (which in surgery is moments)."

"The process for uptake and release of an external substance, such as drugs and contrast agents, are unique in cancerous tissues. As such, we envisaged that an approach combining biophysics-inspired modeling and AI could analyze intraoperative changes in NIR intensities over time in varied tissue, enabling clinically useful lesion classification with high specificity. To translate this knowledge for the first time into an intraoperative surgical decision support tool, a computer vision-AI real time tissue-tracking and categorizing prototype has been developed. As the prototype relies only on the NIR fluorescence data stream, it is usable with commercially available imaging systems."

Also speaking about the publication of this study in the *British Journal of Surgery*, Professor Donal O'Shea, Department of Chemistry at RCSI University of Medicine and Health Sciences said "Targeted agents for cancer imaging currently under trial adhere rigidly to conventional paradigms of fluorescence guided surgery mechanisms, and in the main are administered systemically before surgery, with the operation scheduled for when maximum stable contrast between the tumor and other tissues exist. However this timing is often unpredictable, it can take some days and false positives can occur. Clinical usefulness is further limited by dosing practicalities, scheduling challenges and patient-to-patient and cancer-to-cancer differences. This work instead indicates a

novel pathway and process for immediate, perfect realization of agent information during surgery which would greatly improve efficiency and effectiveness of [cancer](#) care."

More information: Artificial intelligence indocyanine green (ICG) perfusion for colorectal cancer intra-operative tissue classification. *British Journal of Surgery*, znaa004, doi.org/10.1093/bjs/znaa004

Provided by University College Dublin

APA citation: Artificial Intelligence used to enhance decision making during colorectal cancer surgery for first time (2020, December 17) retrieved 18 June 2021 from

<https://medicalxpress.com/news/2020-12-artificial-intelligence-decision-colorectal-cancer.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.